**RESPIRONICS**

# Instructions for Use

## Whisper Swivel II

The Whisper Swivel II is a multi-patient use exhalation port that provides a continuous leak path in the patient circuit when used with positive pressure devices.

**⚠ Caution:** U.S. federal law restricts this device to sale by or on the order of a physician.

**🚫 Note:** This product is not made with natural rubber latex or DEHP.


### Symbols Glossary

Symbol	Title and Meaning	Reference
	Prescription device Caution: U.S. federal law restricts this device by or on the order of a physician.	-
	Caution, consult accompanying documents.	IEC 60878 Symbol 5.1.2, ISO 15223-1
	Batch code Indicates the manufacturer's batch code so that the batch or lot can be identified.	IEC 60878 ISO 7000-2492 Symbol 5.1.5, ISO 15223-1
	Reorder number Indicates the manufacturer's catalogue number so the medical device can be identified.	ISO 7000-2493 Symbol 5.1.6, ISO 15223-1
	Manufacturer Indicates the medical device manufacturer	IEC 60878 ISO 7000-3082 Symbol 5.1.1, ISO 15223-1
	Humidity limitation Indicates the range of humidity to which the medical device can be safely exposed.	IEC 60878 ISO 7000-2620 Symbol 5.3.8, ISO 15223-1
	Temperature limit Indicates the storage temperature limits to which the medical device can be safely exposed.	IEC 60878 ISO 7000-0632 Symbol 5.3.7, ISO 15223-1
	Not Made with Natural Rubber Latex Indicates that this product is not made with natural rubber latex.	-

### ⚠ Warnings:

- The body of the Whisper Swivel II must be inserted into the sleeve until they snap together to ensure a proper fit (Figure 1). Otherwise, in some circumstances, suffocation can occur.
- As with most CPAP devices, at low CPAP pressures, the air flow through the exhalation port may not be enough to clear all the exhaled gas (CO<sub>2</sub>) from the mask. You may rebreathe some of the air that you have exhaled.
- Hand wash prior to first use. Wash and inspect for damage or wear daily. Discard and replace any components as necessary.

- This warning applies to almost all CPAP devices. The Whisper Swivel II should be used with the CPAP system recommended by your healthcare professional or respiratory therapist. As with any exhalation device, do not apply the mask to your face unless the CPAP device is turned on and operating properly. To allow a continuous flow of air out of the mask, the vent holes in the exhalation port must never be blocked. When the CPAP device is turned on and functioning properly, new air from the CPAP device flushes exhaled air out through the exhalation port's vent holes. However, when the CPAP device is not operating, enough fresh air will not be provided through the mask, and exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can, in some circumstances, lead to suffocation.
- Do not tape the Whisper Swivel II parts together. Otherwise, in some circumstances, CO<sub>2</sub> rebreathing can occur.

### Daily Cleaning Instructions for Home Use

1. Disconnect the Whisper Swivel II from the mask and tubing. (Note: To disconnect the tubing, bend at the tubing connection.) Disassemble the Whisper Swivel II by pressing together the two tabs at the top of the swivel separating the two parts.
2. Gently hand wash the Whisper Swivel II parts in a solution of warm water and a mild shampoo or dish detergent.

**⚠ Caution:** Do not use bleach, alcohol, cleaning solutions containing bleach or alcohol, or cleaning solutions containing conditioners or moisturizers.

**⚠ Caution:** Any deviation from the manufacturer's instructions or use of agents other than those listed in this guide may impact the performance or durability of the product.

3. Rinse thoroughly. Allow to air dry.
4. Reassemble the parts by sliding the body into the sleeve until they snap together.
5. Inspect the Whisper Swivel II for any remaining debris or damage. Repeat the cleaning instructions if necessary. Replace if damaged.

### Cleaning and Disinfecting Instructions for Institutional Use

**⊖ Note:** For Institutional Use, the Whisper Swivel II must be disinfected prior to each patient's use. Adequate disinfection can only be performed after the device has been thoroughly cleaned following the pre-processing instructions below.

#### Pre-Processing Instructions

1. Disassemble the Whisper Swivel II by pressing the two tabs at the top of the swivel and separating the two parts. Soak the swivel parts in a commercially available enzymatic cleaner (e.g., Klenzyme™) following the manufacturer's instructions.
2. Manually clean the Whisper Swivel II with a soft bristle brush in a solution of cool tap water and a commercially available anionic detergent (e.g., Manu Klenz™).

**⚠ Caution:** Do not use cleaning solutions containing conditioners or moisturizers.

3. Rinse thoroughly following the manufacturer's instructions. If no rinsing instructions are provided, rinse the swivel thoroughly under cool running tap water. Air dry.

**⊖ Note:** If this procedure does not adequately remove adhering substances, replace the Whisper Swivel II.

4. Reassemble the parts by sliding the body into the sleeve until they snap together.

#### Disinfecting the Whisper Swivel II

1. Disassemble the Whisper Swivel II by pressing the two tabs at the top of the swivel and separating the two parts.
2. Soak the swivel parts in a commercially available 2.4% or 2.5% glutaraldehyde solution (e.g., Cidex™) following the manufacturer's instructions.

- Rinse the swivel parts thoroughly following the manufacturer's instructions. If no instructions are provided, rinse the swivel thoroughly under cool running tap water. Air dry.
- Inspect the swivel parts for signs of wear or deterioration. Replace if damaged.
- Reassemble the parts by sliding the body into the sleeve until they snap together.


**Note:** Testing has shown that the Whisper Swivel II can withstand up to 15 disinfection cycles using a 2.5% glutaraldehyde solution following the manufacturer's instructions.

### Connecting the Whisper Swivel II

- Firmly insert the top end of the Whisper Swivel II into the end of the mask's rotating connector (Figure 2). Slightly twist the Whisper Swivel II to ensure a tight seal. To reduce the amount of air flow blowing onto the user, position the swivel as shown in figure 2.
- Connect the flexible tubing to the rotating end of the Whisper Swivel II.

**Figure 2**


- A = Rotating Connector
- B = Whisper Swivel II
- C = Vent Holes (Exhaled Air)
- D = Flexible Tubing

### Specifications

**Warning:** The technical specifications of the mask are provided for your healthcare professional to determine if it is compatible with your CPAP or bi-level therapy device. If used outside these specifications, or if used with incompatible devices, the mask may be uncomfortable, the seal of the mask may not be effective, optimum therapy may not be achieved, and leak, or variation in the rate of leak, may affect device function. The pressure flow curve shown below is an approximation of expected performance. Exact measurements may vary.

#### Pressure Flow Curve

### Resistance

#### Drop in Pressure at

	50 SLPM	100 SLPM
Drop in Pressure	0.1 cm H <sub>2</sub> O	0.3 cm H <sub>2</sub> O

#### Sound Levels at 10 cm H<sub>2</sub>O

A-weighted Sound Power Level	42.6 dBA
A-weighted Sound Pressure Level at 1m	34.6 dBA

### Disposal

Dispose of in accordance with local regulations.

## Storage Conditions

Temperature: -4° to 140° F (-20° to +60° C)

Relative Humidity: 15% to 95% non-condensing

## Limited Warranty

Respironics, Inc. warrants that the Whisper Swivel II (the "Product") will be free from defects of workmanship and materials for a period of ninety (90) days from the date of purchase (the "Warranty Period"). If the Product fails under normal conditions of use during the Warranty Period and the Product is returned to Philips Respironics within the Warranty Period, as the sole remedy under this warranty, Philips Respironics will replace the Product. This warranty applies only to the original purchaser of the Product and is nontransferable.

### Exclusions

This warranty does not cover damage caused by accident, misuse, abuse, negligence, alteration, failure to use or maintain the Product under conditions of normal use and in accordance with the terms of the product literature, and other defects not related to materials or workmanship.

This warranty does not apply to any Product that may have been repaired or altered by anyone other than Philips Respironics.

### Disclaimers

EXCEPT AS SET FORTH IN THIS LIMITED WARRANTY, PHILIPS RESPIRONICS MAKES NO WARRANTIES, EXPRESSED OR IMPLIED, STATUTORY OR OTHERWISE, REGARDING THE PRODUCT, ITS QUALITY OR PERFORMANCE. PHILIPS RESPIRONICS SPECIFICALLY DISCLAIMS THE IMPLIED WARRANTY OF MERCHANTABILITY AND THE IMPLIED WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL PHILIPS RESPIRONICS BE LIABLE FOR THE COST OF PROCUREMENT OF SUBSTITUTE PRODUCTS OR FOR ANY INDIRECT, SPECIAL, CONSEQUENTIAL, PUNITIVE, EXEMPLARY OR INCIDENTAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, ECONOMIC LOSS, LOST REVENUES, OR OVERHEAD), FOR ANY CAUSE OF ACTION, WHETHER IN CONTRACT OR TORT, AND WHETHER OR NOT PHILIPS RESPIRONICS WAS AWARE OR SHOULD HAVE BEEN AWARE OF THE POSSIBILITY OF THESE DAMAGES. EXCEPT FOR THE OBLIGATIONS UNDER THIS LIMITED WARRANTY, PHILIPS RESPIRONICS WILL HAVE NO OBLIGATION OR LIABILITY FOR ANY OTHER LOSS, DAMAGE OR INJURY RESULTING DIRECTLY OR INDIRECTLY FROM THE PRODUCT. Some jurisdictions do not allow (1) the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you, or (2) disclaimers of implied warranties, so the above disclaimer of implied warranties may not apply to you. This warranty gives you specific legal rights, and you may also have other rights under laws applicable in your specific jurisdiction.

To exercise your rights under this limited warranty, contact your local authorized Respironics, Inc. dealer or Respironics, Inc. at 1001 Murry Ridge Lane, Murrysville, Pennsylvania 15668, USA, or Respironics Deutschland GmbH & Co KG, Gewerbestrasse 17, 82211 Herrsching, Germany.

## Limited Warranty for Australia and New Zealand only

Note: For Australian and New Zealand customers this warranty replaces the warranty contained above.

1. Respironics, Inc. warrants that the Product will be free from defects of workmanship and materials and will perform in accordance with the Product specifications. 2. This warranty is valid for a period of 30 days from the date of purchase from an authorised Philips Respironics dealer. 3. If the Product is found to contain a defect of workmanship

or materials or fails to perform in accordance with the Product specifications during the applicable warranty period, Philips Respironics will repair or replace, at its option, the defective material or part or Product. 4. The customer is responsible for returning the Product to an authorized Philips Respironics dealer, and collecting the product from the authorised Philips Respironics dealer after repair or replacement, at its own cost. Philips Respironics is responsible only for the freight cost of transporting the Product between the authorised Philips Respironics dealer and Philips Respironics. Philips Respironics reserves the right to charge an evaluation and postage fee for any returned Product as to which no problem is found following investigation. 5. This warranty does not cover damage caused by accident, misuse, abuse, alteration, and other defects not related to materials or workmanship. 6. The warranty provided by Philips Respironics herein is not transferrable by the initial Buyer of the Product in the event of any sale or transfer of Products purchased by the Buyer from an authorised Philips Respironics dealer. 7. To the extent permitted by law, where the Buyer has the benefit of an implied guarantee under the Australian Consumer Law, but the Product is not of a kind ordinarily acquired for personal, domestic or household use or consumption Philips Respironics' liability will be limited, at the option of Philips Respironics, to the replacement or repair of the Product or the supply of an equivalent Product. 8. To exercise your rights under this warranty, contact your local authorised Philips Respironics dealer. A list of all authorised dealers is available at the following link: [http://www.healthcare.philips.com/au\\_en/homehealth/distributors\\_index.wpd](http://www.healthcare.philips.com/au_en/homehealth/distributors_index.wpd). Alternatively, you can make a claim under this warranty by contacting Philips Respironics directly at: Philips Electronics Australia Limited, 65 Epping Road, North Ryde NSW 2113, Australia. Tel: 1300 766 488, Email: [prcontact@philips.com](mailto:prcontact@philips.com). 9. The following statement is provided to a Buyer who is a "consumer" under the Australian Consumer Law: Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the good repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. 10. The following statement is provided to a Buyer who is a "consumer" under the Consumer Guarantees Act 1993, New Zealand: Our goods come with guarantees that cannot be excluded under the Consumer Guarantees Act 1993. This guarantee applies in addition to the conditions and guarantees implied by that legislation.

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© 2017 Koninklijke Philips N.V. Sva prava pridržana.
Respironics Inc.
Respironics Deutschland GmbH & Co. KG