Therapy Mask 3100 NC Therapy Mask 3100 SP

Instructions for Use

1 Safety Information

Intended Use

This mask is intended to provide an interface for application of CPAP or bi-level therapy to patients. The mask is for single patient use in the home and multi-patient use in the hospital/institutional environment. The mask is to be used on patients >7 years old (>40 lbs) for whom CPAP or bi-level therapy has been prescribed.

⚠ Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

Note:

- This mask is not made with natural rubber latex or DEHP.
 - An electronic copy of these instructions can be found at: www.philips.com/IFU.
-

Contraindications

Use of this mask is contraindicated in patients with the following metallic implants: brain shunts, aneurysm clips, splinters in the eyes, and certain neurostimulators used in and around the head and neck.

Warning: Indicates the possibility of danger to the user or operator.

Caution: Indicates the possibility of damage to the device.

⚠ Warnings

- Magnets are used in the mask. The magnetic field strength is 400 mT. Ensure the mask is kept at least 2 inches (5 cm) away from any active medical implant or medical device that can be impacted by the magnetic field (e.g., pacemaker, defibrillators, neurostimulators, cochlear implants, hearing aids) to avoid possible effects from localized magnetic fields.
- This mask is not intended for patients who are dependent on mechanical ventilation for their life support.
- Use of this mask while the system is not on and operating may cause rebreathing of exhaled air. Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocation.
- Discontinue use and contact your healthcare professional if any of the following symptoms occur: skin redness, irritation, discomfort, blurred vision, or drying of the eyes.
- Discontinue use and consult your healthcare professional if any of the following symptoms occur: unusual chest discomfort, shortness of breath, severe headache, eye pain or eye infections.

- Discontinue use and consult your healthcare professional if any of the following symptoms occur: tooth, gum, or jaw soreness. Use of a mask may aggravate an existing dental condition.
- Do not block the exhalation ports.
- Do not overtighten the headgear straps. Watch for signs of overtightening, such as excessive redness, sores or bulging skin around the edges of the mask. Loosen the headgear straps to alleviate symptoms.
- The mask contains small parts which could result in a choking hazard.
- Do not use in or near magnetic resonance imaging (MRI) equipment.


Note: Any serious incident that has occurred in relation to this device should be reported to Philips and the competent authority of the Member State in which the user and/or patient is established.

2 Contents

Figure ①

Notes:

- Your Therapy Mask 3100 NC or Therapy Mask 3100 SP comes with either an under the nose nasal cushion (NC) (A2) or a silicone pillows (SP) cushion (A1). These instructions refer to both styles as “the cushion” unless otherwise noted.
- An exhalation port is built into this mask, so a separate exhalation port is not required.

- Ⓐ1 Silicone pillows cushion
- Ⓐ2 Under the nose nasal cushion
- Ⓑ Headgear bottom strap
- Ⓒ Headgear top strap
- Ⓓ Sliders
- Ⓔ Wing
- Ⓕ Mask frame and tube assembly
- Ⓖ Exhalation port (do not block)

Accessories

The following accessories are available for use with this mask. Contact Philips Respironics or your healthcare provider for additional information.

Part #	Description
1006085	Deluxe Chin Strap
1012911	Premium Chin Strap


3 Before Use

- Hand wash the entire mask (see **Mask Care** section).
- Inspect the entire mask. Discard and replace any damaged or worn parts.

4 Achieving the Right Fit

Cushion Sizing

Five cushion sizes (XS, S, M, L, MW) are available. Place the cushion sizing gauge under your nose (as shown). For the best performance, use the smallest cushion that fits your nose.

1. Hold the gauge horizontally under your nose and select the cushion size based on where the outer edge of your nostrils and tip of your nose contact the gauge. A mirror or another person may assist with the sizing.
2. The cushion opening will be directly under your nostrils. The cushion will hug your nose and leak will be minimal.
3. At no time should your nose be inside the cushion opening. If any part of the cushion overlaps your nose, you are wearing the cushion incorrectly. You may need a different cushion size.

Contact your healthcare provider for additional information.



Putting on the Mask

The right and left sides of the headgear fabric are marked with an “L” and an “R.” Before putting on the mask, check that the “L” and the “R” are upright and on the correct sides. ②

1. Wash your face.

⚠ Caution: Do not use moisturizer/lotion on your hands or face.


2. Holding the cushion against your nose with one hand, pull the headgear straps over your head. ③

Adjusting the Mask

1. Use the sliders to adjust the top and bottom headgear straps. To adjust while wearing the mask, grasp one of the sliders in one hand. Grasp the headgear strap fabric in your other hand. Pull the slider to either the left or right to tighten or loosen. ④

Note: The most common mistake is overtightening. If your mask feels tight, loosen the headgear straps to adjust the fit. The mask should fit loosely and comfortably against your face. ⑤

Using the Mask

Note: Re-seat the mask by pulling it away from your face and gently setting it back into place.

1. Connect the flexible tubing (included with your CPAP or bi-level device) to the swivel.
2. Turn the therapy device on. Lie down. Breathe normally.
3. Assume different sleeping positions. Move around until comfortable. If there are any excessive air leaks, make final adjustments. Some air leaking is normal.

5 Mask Care

⚠ Cautions:

- Any deviation from these instructions may impact the performance of the product.
- Do not use bleach, cleaning solutions containing bleach, or cleaning solutions containing conditioners or moisturizers.
- Do not use unapproved cleaners (e.g., ozone gas or UV light products).

Cleaning Instructions for Non-Fabric Mask Parts

⚠ Caution: Do not use alcohol or cleaning solutions containing alcohol on non-fabric mask parts.

Hand wash non-fabric parts daily.

1. Disassemble the mask (see **Disassembly** section).
2. Inspect the mask parts for damage or wear (cracking, crazing, tears etc.). Discard and replace any components as necessary.
3. Mix 1 teaspoon (5 ml) of liquid detergent per 1 gallon (3.8 liters) of warm potable water (80–90°F/27–32°C).
4. Completely submerge the mask parts in the warm water/detergent mixture.
5. While submerged, use a soft bristle brush to scrub the mask parts for at least 5 minutes, taking care to clean any cavities or crevices.

Note: Ensure there are no air pockets present inside the mask parts while submerged.

6. Remove the mask parts from the water/detergent mixture and rinse with 1.5 gallons (5 liters) of water for at least 1 minute.
7. Air dry mask parts out of direct sunlight. Make sure the entire mask is dry before use.

Cleaning Instructions for Fabric Parts

⚠ Caution: Do not place the fabric parts into the dryer.

Hand wash fabric parts weekly, or as needed.

1. Disassemble the mask (see **Disassembly** section).
2. Inspect the fabric parts for damage or wear. Discard and replace any components as necessary.
3. Mix 1 teaspoon (5 ml) of liquid detergent per 1 gallon (3.8 liters) of warm potable water (80–90°F/27–32°C).
4. Completely submerge the fabric parts in the warm water/detergent mixture and soak for 5 minutes.
5. Agitate the fabric parts in the water/detergent mixture for an additional 5 minutes.
6. Remove the fabric parts from the water/detergent mixture and rinse with 1 gallon (3.8 liters) of warm water for at least 1 minute.
7. Lay the fabric parts flat or line dry. Make sure the entire mask is dry before use.

6 Assembly and Disassembly

⚠ Cautions:

- The tubing is permanently attached to the mask frame. Do not attempt to remove the tubing from the frame.
- Do not pull on the mask tube to separate it from CPAP tubing.








Assembly

1. The inside of the wings and the headgear are marked with “R” and “L” to show the right and left sides. Match the left (“L”) side of the headgear to the left wing. Make sure each “L” is facing up. **⑥**
2. Insert the tip of the left wing into the slot in the left side of the headgear fabric. Feed the wing through the length of the fabric. **⑦**
3. Hook the fabric around the edges of the wing, making sure the wing is fully covered. **⑧** If any part of the wing is visible **⊗**, adjust the fabric.
4. Repeat with the right (“R”) side of the wing and headgear.
5. Align the cushion with the frame. The inside of the cushion is marked with an “R” and an “L” to show the left and right sides. Make sure the “R” and “L” are facing up. The magnets in the cushion will snap into place. **⑨**
6. Confirm that the cushion is firmly attached to the frame. **⑩**

Note: If placed in the frame upside down, the cushion will not connect to the magnets.

Disassembly

1. Grasp the back of the cushion and pull it away from the frame. You will feel a pull as the magnets separate. ⑨
2. Pinch the upper corner of the headgear fabric and peel it off the top of the wing. ⑧
Slide the fabric completely off of the wing ⑦, and then repeat with the other side of the headgear.

7 Health Care Provider & Clinician Information

MRI Safety Information

The Therapy Mask 3100 NC/Therapy Mask 3100 SP is MR Unsafe. Keep it outside the MRI scan room (Zone IV). It represents a projectile hazard.

Useful Life

The useful life of the mask depends on the use conditions and maintenance (cleaning, institutional disinfection, and component replacement) of the device. Inspect the mask parts regularly for damage or wear. Replace components as deemed necessary.

Multi-Patient Use

Clean and disinfect the mask daily and between patients when in use in a clinical setting. For multi-patient use in the hospital/institutional environment, use the Disinfection Guide to reprocess the mask between patients. These instructions can be obtained online at www.philips.com/IFU, or by contacting Philips Respironics Customer Service at 1-800-345-6443 (USA or Canada), at +1-724-387-4000, or go to www.respironics.com to find your local customer service contact information.

Philips Respironics DreamStation Mask Type and System One Resistance Control Settings

Using the mask with a Philips Respironics DreamStation or System One device provides optimal comfort. The provider sets this value on the device.

Silicone pillows cushions	
Cushion size	Setting
XS, S	X3
M, MW, L	X2


Under the nose nasal cushions	
Cushion size	Setting
All	X2

Specifications

The technical specifications of the mask are provided for your healthcare professional to determine if it is compatible with your CPAP or bi-level therapy device.

⚠ Warning: Use outside of these specifications may result in ineffective therapy. The pressure flow curve shown below is an approximation of expected performance. Exact measurements may vary.

Pressure Flow Curve

Resistance

Drop in Pressure (cm H ₂ O)		
Under the nose nasal cushion size	at 50 SLPM	at 100 SLPM
XS	1.7	6.3
S	1.5	6.5
M	1.4	6.1
MW	1.6	6.1
L	1.5	5.3

Drop in Pressure (cm H ₂ O)		
Silicone pillows cushion size	at 50 SLPM	at 100 SLPM
XS	1.9	7.5
S	2.0	7.2
M	1.5	6.1
MW	1.6	6.3
L	1.8	6.8

Deadspace

Cushion size	Under the nose nasal cushions	Silicone pillows cushion
XS	14 mL	11 mL
S	18 mL	12 mL
M	18 mL	12 mL
MW	17 mL	13 mL
L	24 mL	12 mL

Sound Levels

DECLARED DUAL-NUMBER NOISE EMISSION VALUES in accordance with ISO 4871

A-weighted sound pressure level @ 1 m: 20 dB(A) with an uncertainty of 6 dB(A)

A-weighted sound power level: 28 dB(A) with an uncertainty of 6 dB(A)

Values determined according to noise test code given in ISO 17510, using the basic standards ISO 3744 and ISO 4871.

Note: The sum of a measured noise emission value and its associated uncertainty represents an upper boundary of the range of values which is likely to occur in measurements.

Disposal

Dispose of the mask and headgear annually.

Dispose of this device in accordance with local collections and recycling regulations. For more information, visit www.philips.com/recycling.



Storage Conditions

Temperature: -4° to 140° F (-20° to +60° C)

Relative Humidity: 15% to 95%

8 Symbols Glossary

Refer to <http://www.symbols.philips.com> for a description of the symbols used on this device and its packaging.

Symbol	Title and Meaning	Symbol	Title and Meaning
	Unique Device Identifier – Indicates the Unique Device Identifier information.		Batch code – Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Humidity limitation – Indicates the range of humidity to which the medical device can be safely exposed.		Temperature limit – Indicates the temperature limits to which the medical device can be safely exposed.
	Catalogue number – Indicates the manufacturer's catalogue number so that the medical device can be identified.	X2 X3	Philips Respironics System One Resistance Control Value – Indicates the therapy device comfort setting.
 <small>www.philips.com/IFU</small>	Electronic instructions for use Indicates that relevant information for use of the product is available in electronic form		Packaging unit – Indicates the number of pieces in the package.
	Caution – Indicates the need for the user to consult the instructions for use.		Manufacturer – Indicates the medical device manufacturer
	Not Made with Natural Rubber Latex		Date of manufacture – Indicates the date when the medical device was manufactured. Country of manufacture – Indicates the country of manufacture of the product. Note: When applied to the label, "CC" is replaced by the country code.
	Prescription device		
	Authorized representative in the European Community		Consult instructions for use
	Importer – Indicates the entity importing the medical device into the locale.		MR unsafe

9 Limited Warranty

Respironics, Inc., a Philips company ("Philips Respironics") provides this non-transferable, limited warranty for the Therapy Mask 3100 NC and Therapy Mask 3100 SP ("Product") to the customer who originally purchased the Product directly from Philips Respironics.

What this Warranty Covers: Philips Respironics warrants each new Product will be free from defects in materials and workmanship and will perform in accordance with the Product specifications under normal and proper use and maintenance in accordance with applicable instructions, subject to the exclusions below.

How Long does this Warranty Last: For the product, ninety (90) days from the longer of the date of shipment to the purchaser or date of setup by purchaser for the end user. The warranty period for accessories, replacement parts, and disposables including, but not limited to, chin straps is ninety (90) days from the date of shipment to the original purchaser.

What this Warranty does not cover: This warranty does not apply to any software included with the Product as the software warranty is included in the software license. This warranty does not cover damage or injury whether to the Products, personal property, or persons caused by accident, misuse, abuse, Acts of God, water ingress, repair or alteration by anyone other than Philips Respironics or its authorized service center, failure to operate in accordance with the terms of the operating manual and instructions, lack of reasonable care, the discontinuance of a network (e.g. 2G, 3G, etc.) by a carrier (e.g. ATT, Verizon, etc., or other defects not related to material or workmanship. This warranty is not transferable. If Philips Respironics finds that a Product returned for service or the issue raised is not covered under this limited warranty, Philips Respironics may charge an evaluation fee and return shipping.

What Philips Respironics will do: If a Product fails to conform to the warranties set forth above during the applicable warranty period, Philips Respironics will repair or replace the Product or refund the original purchase price, in Philips Respironics sole discretion. Philips Respironics may use new or remanufactured assemblies, components, and parts in repair and new or recertified refurbished devices for replacement. The balance of the original warranty period will apply to any Product or component of a Product repaired or replaced under this warranty.

Warranty Disclaimer; Limitation of Liability: EXCEPT AS SET FORTH IN THIS LIMITED WARRANTY, PHILIPS RESPIRONICS MAKES NO WARRANTIES, EXPRESSED OR IMPLIED, STATUTORY OR OTHERWISE, REGARDING THE PRODUCT OR ITS QUALITY OR PERFORMANCE. PHILIPS RESPIRONICS SPECIFICALLY DISCLAIMS THE IMPLIED WARRANTY OF MERCHANTABILITY AND THE IMPLIED WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL PHILIPS RESPIRONICS MAXIMUM LIABILITY UNDER THESE WARRANTIES EXCEED THE ORIGINAL PURCHASE PRICE OR WILL PHILIPS RESPIRONICS BE LIABLE FOR ANY ECONOMIC LOSS, LOSS OF PROFITS, OVERHEAD, OR SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. Repair, replacement, or return of purchase price by Philips Respironics is the original purchaser's sole and exclusive remedy under this warranty.

This warranty gives you specific legal rights, and you may also have other rights that vary from state to state. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above exclusion and limitations may not apply to you.

How to get warranty support: Patients contact your local authorized Philips Respironics dealer and dealers contact Respironics, Inc. at:

1001 Murry Ridge Lane
Murrysville, Pennsylvania 15668-8550
+1-724-387-4000
Deutschland
Gewerbestrasse 17
82211 Herrsching, Germany
+49 8152 93060

Limited Warranty for Australia and New Zealand only

Note: For Australian and New Zealand customers this warranty replaces the warranty contained above.

1. Respireonics, Inc., a Philips company ("Philips Respireonics"), warrants that the Products shall be free from defects of workmanship and materials and will perform in accordance with the Product specifications. 2. This warranty is valid for a period of ninety (90) days from the date of purchase from an authorized Philips Respireonics dealer. 3 If the Product is found to contain a defect of workmanship or materials or fails to perform in accordance with the Product specifications during the applicable warranty period, Philips Respireonics will repair or replace, at its option, the defective material or part. 4. The customer is responsible for returning the product to an authorised Philips Respireonics dealer, and collecting the product from the authorised Philips Respireonics dealer after repair or replacement, at its own cost. Philips Respireonics is responsible only for the freight cost of transporting the product between the authorised Philips Respireonics dealer and Philips Respireonics. Philips Respireonics reserves the right to charge an evaluation and postage fee for any returned Product where no problem is found following investigation. 5. This warranty does not cover:

- products purchased outside of Australia or New Zealand;
- damage caused by accident, misuse, abuse, alteration, pest infestation, liquid ingress, and other defects not related to materials or workmanship.

6. The warranty provided by Philips Respireonics herein is not transferrable by the Buyer in the event of any sale or transfer of Products purchased by the Buyer from an authorized Philips Respireonics dealer. 7. To exercise your rights under this warranty, contact your local authorised Philips Respireonics dealer. A list of all authorised dealers is available at the following link: <http://www.philips.com.au/healthcare>.

Alternatively, you can make a claim under this warranty by contacting Philips Respireonics directly at: Philips Electronics Australia Limited, 65 Epping Road, North Ryde NSW 2113, Australia. Tel: 1300 766 488, Email: repairs-src@philips-easyconnect.com. 8. The following statement is provided to a Buyer who is a "consumer" under the Australian Consumer Law: Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the good repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. 9. The following statement is provided to a Buyer who is a "consumer" under the Consumer Guarantees Act 1993, New Zealand: Our goods come with guarantees that cannot be excluded under the Consumer Guarantees Act 1993. This guarantee applies in addition to the conditions and guarantees implied by that legislation.

AUSTRALIAN SPONSOR DETAILS:

Philips Electronics Australia Ltd.
65 Epping Road, North Ryde, NSW 2113
Australia.

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© Koninklijke Philips N.V., 2021. All rights reserved.
Resironics Inc.
Resironics Deutschland GmbH & Co. KG