MANUEL D'UTILISATION • USER GUIDE • GEBRAUCHSANLEITUNG • MANUAL DE USO • ISTRUZIONI D'USO • INSTRUÇÕES DE USO
BRUKSANVISNING • BRUGERVEJLEDNING • GEBRUIKSAANWIJZING • INSTRUKCJA OBSŁUGI • ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
ИНСТРУКЦИЯ ДЛЯ ПОЛЬЗОВАТЕЛЯ • UŽIVATELSKÝ MANUÁL • 使用说明 • 使用説明書


Respireo SOFT Child non vented

Masque nasal • Nasal mask • Nasenmaske • Mascarilla nasal • Maschera nasale • Máscara nasal • Näsmask • Nesemaske
Nasal maske • Neusmasker • Maska nosowa • Ρινική μάσκα • Носовая маска • Nosní maska • 鼻面罩 • 鼻マスク

FRANÇAIS.....	p. 7
ENGLISH.....	p. 11
DEUTSCH.....	p. 15
ESPAÑOL.....	p. 19
ITALIANO.....	p. 24
PORTUGUÊS.....	p. 28
SVENSKA.....	p. 32
NORSK.....	p. 36
DANSK.....	p. 39
NEDERLANDS.....	p. 43
POLSKI.....	p. 47
ΕΛΛΗΝΙΚΑ.....	p. 52
РУССКИЙ.....	p. 56
ČEŠTINA.....	p. 60
中文.....	p. 63
日本語.....	p. 66

Respireo SOFT Child non vented



1


2


3

4

1 - USAGE PREVU

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ENGLISH

1 - APPLICATION

Respireo SOFT Child non vented non vented is a non-invasive nasal mask used for treating respiratory failure in children aged 1,5 to 6 years.

Respireo SOFT Child non vented is designed for use by babies with controlled non-invasive positive pressure ventilation (NPPV). Respireo SOFT Child non vented is designed for use by a single patient at home or more than one in the hospital environment.

⚠ 2 - IMPORTANT NOTES

- 2.1 The mask must be used under the supervision of an adult who has read and understands the instructions.
- 2.2 The mask must only be used with a ventilator equipped with an active exhalation valve or a respiratory circuit with an exhalation valve.

- 2.3 The operating pressure range is 4-24 cmH₂O.
- 2.4 If you experience discomfort or irritation, or suffer an allergic reaction to any of the components of the mask, please report this to your doctor or respiratory therapist.
- 2.5 If oxygen is used with the device, the oxygen supply must be turned off when it is not in use to prevent oxygen accumulating inside, which is a fire risk.
- 2.6 Never smoke or light naked flames when using oxygen.
- 2.7 Do not use the mask in the presence of nausea or vomiting.
- 2.8 The mask must be cleaned and disinfected. For further details refer to section CLEANING THE MASK.
- 2.9 Check the mask each time before use. If there is any visible damage, report it to the dealer.
- 2.10 The mask must only be used in accordance with these instructions only. No other instructions are applicable.
- 2.11 Some of the components of the mask are small enough to be swallowed by children and must never be left unattended. An adult must always be present.
- 2.12 To ensure proper operation of the hose kit, we recommend you avoid pulling the latter and handle it carefully.

3 - COMPONENTS OF THE MASK (Fig.A)

(a) Main body: frame (1), cushion (2);

(b) Tube assembly: elbow fitting (3), tube (4), swivel fitting (5), cap (6);

(c) Headgear assembly: headgear (7); hooks (8).

4 - HOW TO ASSEMBLE THE MASK (Fig. B)

- 4.1 - Connect the tube assembly (b) by inserting the sphere on the elbow fitting (3) in the slot in the main body (a), as shown in Fig. B1.

4.2 Insert the three straps at the top of the headgear (c) in the three slots in the main body (a), as shown in Fig. B2. Then fasten the straps.

5 - HOW TO PUT THE MASK ON THE CHILD (Fig.C)

5.1 Place the mask gently on the child's face, making sure the cushion (2) fits snugly round the nose. Pull the headgear over the child's head, as shown in Fig. C1.

5.2 Place the hooks (8) in the frame (1) slits as shown in figure C2. Then tighten the straps.

N.B: If needed, Respireo SOFT Child non vented can be installed by removing the hooks (8) from the headgear (7) and inserting the bottom belts of headgear (7) directly into the bottom slots of the main body(a).

5.3 Position the headgear (7) comfortably on the child's head and adjust all the straps, starting from the lower ones, as shown in Fig. C3.

5.4 Connect the swivel fitting (5) to the respiratory circuit of the ventilation unit as shown in figure C4. The child will be more reassured and will accept the treatment more readily if the mask is placed correctly before starting the ventilation unit. Switch on the ventilation unit.

5.5 For optimal airtightness, and to remove any accidental leakage around the bubble, the headgear may need to be adjusted. To do so, gently tighten or loosen the straps starting with the lower straps. Ensure the mask is attached securely to the child's face without the straps being too tight.

N.B. The lower straps only serve to hold the mask in position, so they must be kept fairly loose.

6 - HOW TO REMOVE THE MASK

To remove the mask, unhook one of the two hooks (8).

Lift the headgear over the child's head and remove the mask from the child's nose.

It is not necessary to loosen or undo the other straps (Velcro) from the headgear.

7 - HOW TO DISSAMBLE THE MASK (Fig. D)

7.1 Unfasten and pull out the straps, as shown in Fig. D1, then detach the headgear (7) from the main body (a).

7.2 Detach the tube assembly (b) by pulling the elbow fitting (3) out of the main body (a), as shown in Fig. D2.

8 - HOW TO CLEAN THE MASK

8.1 – DAILY CLEANING FOR REUSE BY THE SAME PATIENT

It is recommended to clean the main body daily (a) and the headgear (7) and hooks (8) once a week.

First wash your hands thoroughly.

Disassemble all the components (see HOW TO DISASSEMBLE THE MASK) and wash carefully in warm water at about 30°C, using mild soap. Rinse thoroughly under the tap and leave in the air to dry, but not in direct sunlight.

8.2 – REPROCESSING THE MASK FOR RE-USE WITH OTHER PATIENTS

N.B. The mask must always be regenerated before being used by another patient.

N.B. The headgear (c) must be replaced each time.

N.B. In case of thermal disinfection and autoclave sterilization, the tube assembly (b) must be replaced.

The mask can be regenerated up to 12 times, using one of the following methods:

- high-level thermal disinfection;
- high-level chemical disinfection;
- sterilisation in an autoclave at 121°C for no more than 20 minutes;

Full details of the mask regeneration process can be obtained by logging on to www.device.airliquidehealthcare.com or by contacting your dealer.

9 - PRECAUTIONS

Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents or aromatic oils as they would damage the product and reduce its life.

Prolonged exposure to direct sunlight or heat may damage the mask.

If any damage is noticed (cracks, tears, etc.), the damaged part must be replaced.

Check carefully to make sure there are no soap deposits or foreign bodies inside the mask.

10 - PRESSURE/OXYGEN PORT

If additional oxygen or a pressure reading is required, use the port in the elbow fitting (3).

Remove the cap (6) and insert the oxygen/pressure tube in the port.

11 - TROUBLE-SHOOTING CHART

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
<ul style="list-style-type: none"> • No air or insufficient air comes out of the mask 	<ul style="list-style-type: none"> • Flow generator off • Flow generator regulated wrongly • Patient tube not connected 	<ul style="list-style-type: none"> • Switch on the flow generator. • Regulate the flow generator. • Connect the patient tube to the flow generator and mask.
<ul style="list-style-type: none"> • The mask is noisy 	<ul style="list-style-type: none"> • Mask assembled incorrectly 	<ul style="list-style-type: none"> • Take the mask apart and reassemble it, as explained in section HOW TO ASSEMBLE THE MASK.
<ul style="list-style-type: none"> • The mask does not fit the face snugly 	<ul style="list-style-type: none"> • Mask not positioned correctly on the patient's face • Headgear straps too tight 	<ul style="list-style-type: none"> • Reposition the mask, as explained in section HOW TO PUT THE MASK ON THE CHILD. • Unfasten the lower straps and refasten more loosely.

12 - DISPOSAL INSTRUCTIONS

The components of the mask are made of the following materials and must be disposed of in accordance with the applicable laws.

COMPONENT	MATERIAL
Main body (a)	Silicone + Polycarbonate
Tube assembly (b)	Polycarbonate + Metallocene + Acetalic + Silicone
Headgear assembly (c)	Nylon + polyurethane

13 - TECHNICAL DATA

Normative references	CE marking in accordance with European Directive 93/42/EEC, class IIa, as amended, EN ISO 14971, ISO 5356-1, ISO 10993-1 and EN ISO 17510-2, EN ISO 17664
Dead space volume inside the mask	44 ml
Operating pressure	4-24 cmH ₂ O
Temperature range	5-40°C
Max. overall dimensions	89 mm x 122 mm x 222 mm

14 - WARRANTY

Air Liquide Medical Systems S.r.l. warrants the product free from material and manufacturing defects for a period of 120 days from the date of purchase, provided the operating conditions and instructions are complied with.

If the product is found to be faulty during normal use, it will be repaired or replaced at the manufacturer's discretion.

The manufacturer accepts liability for the safety, reliability and performance of the device only if it is used for the intended application and in accordance with the operating instructions.

The manufacturer accepts liability for the safety, reliability and performance of the device only if repairs are carried out by the manufacturer or an authorised workshop.

The original packaging must be retained for the entire warranty period.

Claims under the warranty should be submitted to your local dealer, who can also provide further details of your warranty rights.

COPYRIGHT

The information contained in this booklet must not be used for any purpose other than the reason for which it is provided. This booklet is the property of Air Liquide Medical Systems S.r.l. and no part of it may be reproduced without the manufacturer's written permission. All rights reserved.

TECHNICAL UPDATES

Air Liquide Medical Systems S.r.l. periodically reviews and modifies all its medical devices in order to improve their performance, safety and reliability. The instruction booklets are updated to make sure they remain consistent with the features of the devices placed on the market. If the booklet accompanying the device gets damaged or mislaid, a replacement copy can be obtained from the manufacturer by quoting the data shown on the rating plate.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.