MANUEL D'UTILISATION • USER GUIDE • GEBRAUCHSANLEITUNG • MANUAL DE USO • ISTRUZIONI D'USO • INSTRUÇÕES DE USO • BRUKSANVISNING
GEBRUIKSAANWIJZING • INSTRUKCJA OBSŁUGI • ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ • ИНСТРУКЦИЯ ДЛЯ ПОЛЬЗОВАТЕЛЯ • UŽIVATELSKÝ MANUÁL • 使用说明 • 使用説明書

Respireo Primo P


MASQUE NARINAIRE avec fuites calibrées • NASAL PILLOW MASK vented • NASENPOLSTER-MASKE mit kalibrierten Luftauslassöffnungen
• MASCARILLA CON ALMOHADILLAS NASALES con orificios de exhalación calibrados • MASCHERA A CUSCINETTI NASALI con fori di
esalazione calibrati • MÁSCARA COM ALMOFADAS NASAIS com furos de exalação calibrados • NÁSMASK MED MJUKDEL med utandningshål
• NEUSMASKER met gekalibreerde ventilatieopeningen • MASKA NOSOWA Z PODUSZKAMI ORAZ z kalibrowanymi otworami wentylacyjnymi •
ΡΙΝΙΚΗ ΜΑΣΚΑ με ρυθμισμένες οπές εξαερισμού • МАСКА С НОСОВЫМИ НАСАДКАМИ с калиброванными выдыхательными отверстиями
• MASKA S NOSNÍMI VLOŽKAMI s kalibrovanými vydechovacími otvory • 鼻枕面罩 带有标准通气口 • 鼻小クッション付きマスク 通気穴付き


FRANÇAIS	p. 6
ENGLISH	p. 9
DEUTSCH	p. 12
ESPAÑOL	p. 16
ITALIANO	p. 20
PORTUGUÊS	p. 24
SVENSKA	p. 27
NEDERLANDS	p. 30
POLSKI	p. 34
ΕΛΛΗΝΙΚΑ	p. 38
РУССКИЙ	p. 41
ČEŠTINA	p. 45
中文	p. 48
日本語	p. 51

Respireo Primo P


A


1


2

3


4

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
C


②


D

①

②

1 - UTILISATION

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1 - APPLICATION

Respireo Primo P is a nasal pillow mask with calibrated vent ports (Vented) for use by a single patient.

Respireo Primo P is designed for use by patients who have been prescribed continuous positive airway pressure (CPAP) or Bi-level positive airway pressure (BiPAP) treatment at home or a hospital.

Respireo Primo P is suitable for adult patients (>30 kg).

⚠ 2 - IMPORTANT NOTES

- 2.1 The mask must only be used with CPAP or Bi-level flow generators recommended by a doctor or respiratory therapist.
- 2.2 The operating pressure range is 4 - 40 cmH₂O.
- 2.3 The masks must only be used with the appliance CPAP or Bi-level on and operating at a pressure of at least 4 cmH₂O. With low pressures using a CPAP appliance, the flow through the vent ports may be insufficient to eliminate all the CO₂ from the mask and the tube, which could lead to partial re-inhalation.
- 2.4 The vents in the mask must never be obstructed.
- 2.5 If you experience discomfort or irritation, or suffer an allergic reaction to any of the components of the mask, please report this to your doctor or respiratory therapist.
- 2.6 If oxygen is used, the supply must be turned off when the appliance is not in operation, to prevent oxygen accumulating inside, which is a fire risk.
- 2.7 Do not smoke or light naked flames when using oxygen.
- 2.8 Do not use the mask in the presence of nausea or vomiting.
- 2.9 The mask must be cleaned and/or disinfected. For further details refer to section 7 - CLEANING THE MASK.
- 2.10 Check the mask before use. If there is any visible damage, report it to the dealer.
- 2.11 The mask must only be used in accordance with these instructions. No other instructions are applicable.
- 2.12 Some of the components of the mask are small enough to be swallowed by children and must never be left unattended. An adult must always be present.

3 - COMPONENTS OF THE MASK (Fig. A)

a) nasal pillow; b) main body; c) tube assembly; d) tube strap; e) rear buckle; f) top buckle.

REMARK: The pack contains three nasal pillows sizes S, M and L.

4 - HOW TO PUT ON THE MASK (Fig. B)

The mask comes fully assembled (size M) and ready for use.

Proceed as follows to put on the mask.

Insert the tube assembly (c) in the tube strap (d).

Hold the main body (b) firmly in one hand, keeping the sidepieces upwards. With the other hand, extend the stretch part of the main body with the other hand (Fig. B1).

Place the main body (b) over your head and insert the nasal pillows (a) into your nostrils (Fig. B2).

Adjust the tension of the rear strap (Fig. B3).

Adjust the tension of the top strap (Fig. B4).

Connect the mask to a CPAP or Bi-level appliance.

5 - HOW TO REMOVE THE MASK

To remove the mask, pull the nasal pillows (a) out of the nose and lift the main body (b) off your head. You do not need to loosen the buckles.

6 - HOW TO DISASSEMBLE THE MASK (Fig. C)

The mask comes fully assembled (size M) and ready for use.

When it needs to be cleaned, proceed as follows to disassemble it.

Remove the elbow from the main body (b) to disconnect the tube assembly (c) from the mask (Fig. C1).

Remove the nasal pillows (a) from the main body (b) (Fig. C2).

7 - CLEANING OF THE MASK

The mask must be washed before each use, so it needs to be disassembled (see 6 – HOW TO DISASSEMBLE THE MASK).

DAILY CLEANING

Wash the nasal pillows (a), the tube assembly (c) and the plastic part of the main body (b) in lukewarm water (about 30 °C) with mild soap. Rinse thoroughly in water and allow to dry naturally.

WEEKLY CLEANING

Wash the nasal pillows (a), the tube assembly (c) and the main body (b), including the fabric part, as instructed under daily cleaning.

The whole mask can be cleaned in a washing machine at 30 °C, but without spin-drying (delicate wash) up to 20 times.

8 - PRECAUTIONS

Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents or aromatic oils. Such solutions would damage the product and reduce its life.

Prolonged exposure to direct sunlight or heat may cause damage to the mask.

If any visible damage is noted (cracks, tears, etc.), the damaged part must be replaced.

Check carefully to make sure there are no soap deposits or foreign bodies inside the mask.

Check regularly that the vents are not obstructed. Clean with a soft brush if necessary.

9 - HOW TO REASSEMBLE THE MASK (Fig. D)

Insert the nasal pillows (a) in the slot in the main body (b). Make sure the two reference points are aligned, then push the nasal pillows as shown in the diagram (Fig. D1).

Insert the elbow of the tube assembly (c), in the slot in main body (b) (Fig. D2).

If you wish, you can secure the tube (c) in the tube strap (d).

10 - TROUBLE-SHOOTING CHART

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
• No air or insufficient air comes out of the mask	<ul style="list-style-type: none"> • Flow generator off • Flow generator regulated wrongly • Patient circuit tube not inserted 	<ul style="list-style-type: none"> • Switch on the flow generator • Regulate the flow generator • Connect the tube to the flow generator and mask
• The mask is noisy	<ul style="list-style-type: none"> • Nasal pillows mounted incorrectly and air leak 	<ul style="list-style-type: none"> • Mount the cushions correctly as instructed
• The nasal pillows do not fit into the nasal cavity	<ul style="list-style-type: none"> • You are using the wrong size 	<ul style="list-style-type: none"> • Replace with a more suitable size (sizes available S - M - L)
• Exhaled air does not flow out and difficulty inhaling	<ul style="list-style-type: none"> • Vent ports obstructed 	<ul style="list-style-type: none"> • Remove the obstruction


11 - DISPOSAL INSTRUCTIONS

The components of the mask are made of the following materials and must be disposed of in accordance with the applicable laws.

COMPONENT	MATERIAL
Nasal pillows	Silicone
Main body	Nylon
Elbow	Polycarbonate
Tube	Metalocene
Connector	Acetalic
Swivel fitting	Polycarbonate
Headgear	Nylon/Polyurethane
Tube strap	Nylon
Rear buckle	Nylon
Top buckle	Nylon

12 - TECHNICAL DATA

MINIMUM CONTINUOUS FLOW

REMARK: Due to manufacturing differences, the exhalation flow rate may differ from the values shown in the table ($\pm 10\%$).

Normative references	CE marking in accordance with European Directive 93/42/EEC, class IIa and subsequent amendments. EN ISO 14971, ISO 5356-1, ISO 10993-1, EN ISO 17510-2
Resistance	at 50 lpm: 0,6 cmH ₂ O - at 100 lpm: 2 cmH ₂ O
A-weighted sound pressure level according to ISO 3744	11,5 dB(A) \pm 1,5 dB(A)
A-weighted sound power level according to ISO 3744	18,2 dB(A) \pm 1,5 dB(A)
Operating pressure	4 – 40 cmH ₂ O
Temperature range	5 – 40 °C
Dead space (mask and tube assembly)	81 ml

13 - LIMITED WARRANTY

This product by Air Liquide Medical Systems S.r.l. is guaranteed against material and manufacturing defects for a period of 120 days as from the date of purchase, provided it is used in accordance with the operating instructions and for the intended purpose only.

If the product fails under conditions of normal use, Air Liquide Medical Systems S.r.l. will repair or replace, at its option the defective product or any of its components.

The manufacturer is liable for the safety, reliability and performance of the apparatus, provided that it is used in accordance with the instructions and for the intended purpose only, and that any repairs are made by the manufacturer or an approved service centre.

Throughout the warranty period keep the original packing.

Warranty claims on effective product must be made by the consumer contacting the dealer.

For further information on your warranty rights contact your dealer.

COPYRIGHT

The information contained in this booklet may not be used for anything other than the intended purpose. This booklet is the property of Air Liquide Medical Systems S.r.l. and no part of it may be reproduced without the manufacturer's written permission. All rights reserved.

TECHNICAL UPDATES

Air Liquide Medical Systems S.r.l. periodically reviews and modifies all its medical devices in order to improve their performance, safety and reliability. The instruction booklets are updated to ensure that they continue to be consistent with the characteristics of the devices marketed. If the booklet accompanying the apparatus is mislaid, a replacement copy can be obtained from the Manufacturer by quoting the data shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.