1. MANUEL D'UTILISATION 2. USER GUIDE 3. GEBRAUCHSANLEITUNG 4. MANUAL DE USO 5. ISTRUZIONI D'USO 6. INSTRUÇÕES DE USO
7. BRUKSANVISNING 8. GEBRUIKSAANWIJZING 9. INSTRUKCJA OBSŁUGI 10. ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ 11. ИНСТРУКЦИЯ ДЛЯ ПОЛЬЗОВАТЕЛЯ
12. UŽIVATELSKÝ MANUÁL 13. 用戶手冊 14. 使用说明 15. 使用説明書

Respireo **SOFT**


Facial Non Vented

1. MASQUE BUCCO-NASAL 2. FACIAL MASK 3. MUND-NASEN-MASKE 4. MASCARILLA ORO-NASAL 5. MASCHERA ORO-NASALE NON VENTED
6. MÁSCARA ORO-NASAL 7. HELMmask 8. MOND-NEUSMASKER 9. MASKA USTNO-NOSOWA 10. ΠΙΝΙΚΗ ΣΤΟΜΑΤΙΚΗ ΜΑΣΚΑ
11. ПОТО-НОСОВАЯ МАСКА 12. NOSO-ÚSTNÍ MASKA 13. 口鼻面膜 14. 口鼻面罩 15. 鼻口腔マスク




1. FRANÇAIS.....	p. 17
2. ENGLISH.....	p. 19
3. DEUTSCH.....	p. 22
4. ESPAÑOL.....	p. 24
5. ITALIANO	p. 27
6. PORTUGUÊS.....	p. 29
7. SVENSKA	p. 32
8. NEDERLANDS.....	p. 34
9. POLSKI.....	p. 36
10. ΕΛΛΗΝΙΚΑ.....	p. 39
11. РУССКИЙ.....	p. 41
12. ČEŠTINA.....	p. 44
13. 中國傳統.....	p. 46
14. 中文	p. 48
15. 日本語.....	p. 50

A


B

1. COMMENT PORTER LE MASQUE 2. HOW TO WEAR THE MASK 3. KORREKTES ANLEGEN DER MASKE 4. CÓMO UTILIZAR LA MÁSCARA
5. COME INDOSSARE LA MASCHERA 6. COME VESTIR A MÁSCARA 7. TA PÅ MASKEN SÅ HÅR 8. HET MASKER OPZETTEN 9. JAK ZAŁOŻYĆ
MASKĘ 10. ΠΩΣ ΦΟΡΑΜΕ ΤΗ ΜΑΣΚΑ 11. ΚΑΚ НАДЕВАТЬ МАСКУ 12. JAK MASKU NASADIT 13. 如何佩戴鼻罩 14. 如何佩戴鼻罩 15. マスクの着用方法

1**2****3**

B


4

5


6





1. COMMENT AJUSTER LE MASQUE 2. HOW TO ADJUST THE MASK 3. KORREKTES EINSTELLEN DER MASKE 4. CÓMO REGULAR LA MÁSCARA 5. COME REGOLARE LA MASCHERA 6. COME REGULAR A MÁSCARA 7. JUSTERA MASKEN SÅ HÄR 8. HET MASKER REGELEN 9. JAK WYREGULOWAĆ MASKE 10. ΠΩΣ ΡΥΘΜΙΖΟΥΜΕ ΤΗ ΜΑΣΚΑ 11. ΚΑΚ ΡΕΓΥΛΙΡΟΒΑΤΉ ΜΑΣΚΥ 12. JAK MASKU UPRAVIT REGELEN 13. 如何調整鼻罩 14. 如何調整鼻罩 15. マスクの調整方法


1. Remarque : ne serrez pas trop les lanières 2. Note: do not over-tighten the straps 3. Hinweis: Ziehen Sie die Gurte nicht zu fest an 4. Nota: no ajustar excesivamente las correas 5. Nota: non stringere eccessivamente le cinghie 6. Nota: não apertar excessivamente as correias 7. Anmärkning: dra inte åt remmarna för hårt 8. Opmerking: de riempjes niet te veel aanspannen 9. Uwaga: nie należy nadmierne dociskać pasków 10. Σημείωση: μην σφίγγετε υπερβολικά τις ζώνες 11. Примечание: не затягивайте чрезмерно ремни 12. Poznámka: popruhy příliš neutahujte 13. 備註: 勿過度收緊繫帶 14. 注: 不要將皮帶拉得過緊 15. 備考: ベルトの締め付け過ぎに注意してください。

D



1. COMMENT SE PRÉPARER À LA THÉRAPIE 2. HOW TO PREPARE FOR TREATMENT 3. VORBEREITUNG AUF DIE THERAPIE 4. CÓMO PREPARARSE PARA LA TERAPIA 5. COME PREPARARSI ALLA TERAPIA 6. COMO PREPARAR-SE À TERAPIA 7. FÖRBERED DIG FÖR TERAPIN SÅ HÄR 8. ZICH VOORBEREIDEN OP DE THERAPIE 9. JAK PRZYKOTOWAĆ SIĘ DO TERAPII 10. ΠΡΟΕΤΟΙΜΑΣΙΕΣ ΘΕΡΑΠΕΙΑΣ 11. КАК ПОДГОТОВИТЬСЯ К ТЕРАПИИ 12. JAK SE PŘIPRAVIT NA TERAPII 13. 如何準備接受治療 14. 如何准备治疗 15. 治療への準備方法

4


5

**1. Remarques 2. Notes 3. Hinweise 4. Notas 5. Note 6. Notas 7. Anmärkning 8. Opmerkingen 9. Uwagi 10. Σημειώσεις
11. Примечания 12. Poznámky 13. 注意 14. 注意 15. 注意事項**

- ▶ **1.** En cas de plis sur la bulle (a), répétez l'opération **2.** If the cushion (a) is wrinkled, repeat the operation **3.** Wenn das Kissen (a) Falten aufweist, wiederholen Sie den Vorgang **4.** Si el cojín (a) presenta arrugas, repita la operación **5.** Se il cuscino (a) presenta grinze, ripetere l'operazione **6.** Se a almofada (a) estiver enrugada, repetir a operação **7.** Om dynan (a) uppvisar skrynklor, ska du upprepa ingreppet **8.** Als plooiën op het kussen (a) aanwezig zijn, de handeling herhalen **9.** Jeśli poduszka (a) marszczy się, powtórzyć czynność **10.** Αν το μαξιλαράκι (a) παρουσιάζει ζάρες, επαναλάβετε τη διαδικασία **11.** Если на подушке (a) имеются складки (a), повторите операцию **12.** Pokud polštářek (a) vykazuje záhyby, operaci opakujte **13.** 如果墊子 (a) 上有褶皺，請重複操作 **14.** 如衬墊(a) 有皺紋，請重複操作 **15.** クッション(a) にしわが寄っている場合は、この作業を繰り返してください
- ▶ **1.** En cas de fuite, resserrez légèrement les lanières du harnais (d) **2.** If there are leaks, tighten the headgear straps (d) **3.** Wenn es zu Undichtigkeiten kommt, ziehen Sie die Gurte des Kopfteils leicht an (d) **4.** En caso de pérdidas, ajuste ligeramente las correas de la cofia (d) **5.** In caso di perdite, stringere leggermente le cinghie della cuffia (d) **6.** Em caso de vazamentos, apertar ligeiramente as correias do protetor de cabeça (d) **7.** Vid läckage, dra åt hättans remmar något (d) **8.** Bij verlies de riempjes van de kap (d) enigszins aanspannen **9.** W przypadku wycieków, lekko docisnąć pasy uprząży (d) **10.** Σε περίπτωση διαρροών, σφίξτε ελαφρώς τις ζώνες του ακουστικού (d) **11.** в случае утечек слегка затяните ремни пращи (d) **12.** V případě úniků jemně přitlačte popruhy krytu (d) **13.** 如果出現洩漏，輕輕地收緊頭套 (d) 的繫帶 **14.** 如有漏气现象，请轻轻拉紧鼻罩的皮带(d) **15.** 漏れが見つかった場合は、ヘッドギア(d)のベルトを軽く締め付けてください

1. COMMENT ENLEVER ET DÉMONTÉ LE MASQUE 2. HOW TO REMOVE AND DISASSEMBLE THE MASK 3. KORREKTES ABNEHMEN UND DEMONTIEREN DER MASKE 4. CÓMO EXTRAER Y DESMONTAR LA MÁSCARA 5. COME RIMUOVERE E SMONTARE LA MASCHERA 6. COMO REMOVER E DESMONTAR A MÁSCARA 7. TA AV OCH MONTERA NED MASKEN SÅ HÅR 8. HET MASKER VERWIJDEREN EN DEMONTEREN 9. JAK ZDJAĆ I ZDEMONTOWAĆ MASKĘ 10. ΠΩΣ ΑΦΑΙΡΟΥΜΕ ΚΑΙ ΑΠΟΣΥΝΑΡΜΟΛΟΓΟΥΜΕ ΤΗ ΜΑΣΚΑ 11. ΚΑΚ ΣΧΙΜΑΤΉ Ι ΡΑΖΒΙΡΑΤΉ ΜΑΣΚΥ 12. JAK MASKU ODSTRANIT A DEMONTOVA 13. 如何取下和拆卸鼻罩 14. 如何移除和拆卸鼻罩 15. マスクの取り外しと解体方法

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1 - INTENDED USE

Respireo SOFT Facial Non Vented is a reusable oro-nasal mask without calibrated exhalation orifices designed to be used by a single patient at home or by multiple patients in a hospital setting.

Respireo SOFT Facial Non Vented is intended to be used by adult patients (weighing over 30 kg), who have been prescribed non-invasive treatment in controlled ventilation therapy (NIPPV).

⚠ WARNINGS

- The mask should be used exclusively with ventilators equipped with an active exhalation valve and/or a breathing circuit with an exhalation valve.
- In case of discomfort, irritation or allergic reactions to any component of the mask, consult your doctor or respiratory physiotherapist.
- In case of additional oxygen administration, smoking or using open flames is prohibited.
- When using oxygen and the positive pressure equipment is not in operation, turn off the oxygen dispenser to prevent accumulation resulting in a risk of fire.
- Do not use the mask if there are signs of vomiting or nausea.
- Keep away from light.
- The mask should be cleaned and/or disinfected (see illustrations in section F).
- Before use, check the integrity of the mask. If damaged during transport, notify your dealer.
- In case of visible deterioration (cracks, tears, etc.), it is necessary to eliminate and replace the damaged mask component.
- Do not leave the components of the mask unattended, some of these may be swallowed by children.
- It is recommended to wash the mask before the first use.
- Do not disassemble the blue elbow fitting (c) from the mask body (b) under any circumstances, even when cleaning the mask.

2 - MASK COMPONENTS (Section A)

a) cushion; b) mask body; c) blue elbow fitting; d) swivel fitting; e) headgear; f) buckle.

3 - MASK USE

To use the mask correctly, see the illustrations in Sections B — C — D — E — G.

4 - REGENERATION OF THE MASK FOR REUSE ON DIFFERENT PATIENTS

Notes:

- The mask should always be regenerated before being used on a new patient.

- The headgear (e) must be replaced.

The mask can be regenerated for up to 20 cycles by one of the following processes:

- thermal disinfection at a high level;
- chemical disinfection at a high level;
- autoclave sterilization.

Detailed mask regeneration instructions can be found on the Internet at www.device.airliquidehealthcare.com or from your retailer.

5 - TROUBLESHOOTING

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
Difficulty breathing.	Patient circuit not inserted correctly.	Connect the patient circuit
Noisy mask.	<ul style="list-style-type: none"> ● The mask is assembled incorrectly; ● The mask is not positioned correctly. 	<ul style="list-style-type: none"> ● Disassemble the mask and then reassemble it as described in the user manual. ● Reposition the mask on your face and adjust the headgear as described in the user manual.

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
Air leaks between the cushion and the face.	<ul style="list-style-type: none"> ● The mask is not properly positioned on the face. ● The straps of the headgear are not properly adjusted (too tight or too loose). ● The size of the mask may not be correct. 	<ul style="list-style-type: none"> ● Reposition the mask and adjust the headgear as described in the user manual. ● Contact your doctor or the technician to find the correct size

6 - DISPOSAL

The mask components are made of the materials indicated below. They do not contain hazardous substances and can therefore be disposed of with ordinary household waste.

COMPONENT	MATERIAL
Cushion (a)	Silicone + Polycarbonate
Mask body (b)	Polycarbonate
Blue elbow fitting (c)	Polycarbonate
Headgear (d)	Nylon + Polyurethane
Buckle (e)	Nylon

7 - TECHNICAL DATA

Standards	CE marking in accordance with European Directive 93/42/EEC, Class IIa and subsequent implementation modifications - UNI EN ISO 14971, UNI EN ISO 5356-1, UNI EN ISO 10993-1, UNI EN ISO 17664
Treatment pressure	4 cmH ₂ O – 40 cmH ₂ O

Dead space (size L)	290 ml
Overall dimensions maximum (without headgear)	174 mm (A) x 116 mm (L) x 115 mm (P)
Weight (size L)	113 g

8 - WARRANTY CERTIFICATE

Air Liquide Medical Systems S.r.l. guarantees its product against material or manufacturing defects for 180 days from the date of purchase, provided that the conditions of use in the instructions for use are observed.

If the product is faulty under normal use, Air Liquide Medical Systems S.r.l. will, at its discretion, repair or replace the faulty product or its components.

The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the operating instructions and its intended use. The manufacturer shall be deemed responsible for the safety, reliability and performance of the device if repairs are carried out at the manufacturer or at a laboratory authorized by the manufacturer.

Keep the original casing for the entire warranty period. The warranty in case of product defects can be applied by contacting your dealer. For more information on the rights under the warranty, please contact your retailer.

COPYRIGHT

All the information in this manual may not be used for purposes other than the original. This manual is the property of Air Liquide Medical Systems S.r.l. and it may not be reproduced, in whole or in part, without written permission from the company. All rights are reserved.

TECHNICAL CHARACTERISTICS UPDATES

For a continuous improvement of performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l.

are regularly subjected to review and change.

Instruction manuals are therefore modified to ensure their constant consistency with the characteristics of devices placed on the market. If the instruction manual accompanying this device is lost, contact the manufacturer to get a copy of the version corresponding to the device by quoting the references shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.