
Respireo NINA

Masque nasal avec fuites calibrées • Nasal mask vented • Nasenmaske mit kalibrierten luftauslassöffnungen • Mascarilla nasal con orificios de exhalación calibrados • Maschera nasale con fori di esalazione calibrati • Máscara nasal com furos de exalação calibrados
Näsmask med kalibrerade utandningshål • Neusmasker met gekalibreerde ventilatieopeningen • Maska nosowa z kalibrowanymi otworami wentylacyjnymi • Ρινική μάσκα με ρυθμισμένες οπές εξαερισμού • Носовая маска с калиброванными выдыхательными отверстиями • Nosní maska s kalibrovanými vydechovacími otvory • 带有标准通气口的鼻面罩 • 鼻マスク 通気穴付き


FRANÇAIS.....	p. 11
ENGLISH.....	p. 14
DEUTSCH.....	p. 17
ESPAÑOL.....	p. 21
ITALIANO.....	p. 25
PORTUGUÊS.....	p. 28
SVENSKA.....	p. 32
NEDERLANDS.....	p. 35
POLSKI.....	p. 39
ΕΛΛΗΝΙΚΑ.....	p. 42
РУССКИЙ.....	p. 46
ČEŠTINA.....	p. 50
中文.....	p. 54
日本語.....	p. 57

A

B


①

②


③


B


④

⑤


⑥


C


1


2


3


4


5





6

D

E**1****2****3****4**

1 - DESTINATION D'USAGE

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1 – INTENDED USE

Respireo NINA is a reusable nasal mask with calibrated exhalation orifices (Vented) designed to be used by a single patient at home or by multiple patients in a hospital setting.

Respireo NINA is intended to be used by adult patients (weighing over 30 kg), who have been prescribed non-invasive positive pressure ventilation therapy (NIPPV), e.g. continuous positive airway pressure therapy (CPAP) or Bilevel (positive airway pressure therapy - BiPAP).

⚠ WARNINGS

- The mask should be used exclusively with positive pressure equipment recommended by your doctor or respiratory therapist.
- The mask should only be used with the positive pressure equipment in operation. Exhaled air evacuation through the exhalation orifices is guaranteed if the positive pressure apparatus supplies a pressure of no less than 4 cmH₂O. At a pressure lower than 4 cmH₂O, the amount of air coming out from the exhalation orifices will not be enough and partial re-inhalation may occur.
- Never cover the exhalation orifices. Before using the mask, always check that the orifices are not clogged, even partially, as this may lead to partial re-inhalation.
- In case of discomfort, irritation or allergic reactions to any component of the mask, consult your doctor or respiratory physiotherapist.
- In case of additional oxygen administration, smoking or using open flames is prohibited.
- When using oxygen and the positive pressure equipment is not in operation, turn off the oxygen dispenser to prevent accumulation resulting in a risk of fire.
- Do not use the mask if there are signs of vomiting or nausea.
- Keep away from light.
- The mask should be cleaned and/or disinfected.
- Before use, check the integrity of the mask. If damaged during

transport, notify your dealer.

- In case of visible deterioration (cracks, tears, etc.), it is necessary to eliminate and replace the damaged mask component.
- Do not leave the components of the mask unattended, some of these may be swallowed by children.
- It is recommended to wash the mask before the first use.
- Contact your dealer as a healthcare professional in case of functional problems, size problems or for any doubt concerning the medical device application.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

2 - MASK COMPONENTS (Section A)

a) cushion; b) mask body; c) elbow fitting; d) swivel fitting; e) headgear; f) buckle.

3 - MASK FITTING (Section B)

- For the first installation, first of all remove cushion (a) from its transparent bag and then attach it to the mask frame (b) by coupling initially the upper part of the cushion to the superior tooth of the frame, then put the fingers and push in order to attach also the inferior part of the cushion to the frame (B.1).
- Keep the mask in one hand and with the other hand pull the headgear (e) over and behind the head (B.2).
- Connect the two lower straps of the headgear (e) to the hooks. Ensure the lower straps pass below the ears (B.3).
- Gently tighten the straps to ensure stability of the mask according to the sequence of B.4.
- Connect the swivel (d) to the device tubing and then switch on the device at the prescribed pressure (B.5).
- Lie down and, if necessary, readjust the forehead straps (B.6).

4 - MASK USE

To disassemble, wash and assemble the mask correctly, see the illustrations in Sections C - D - E.

PRECAUTIONS

- Do not remove the mask without first disconnecting the lower headgear straps (e).
- Wash the headgear (e) and buckle (f) every week.
- Be very careful to remove any dirt from the exhalation holes in the elbow fitting (c).
- Dry the mask components in the air away from direct sunlight.
- Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents, hydrating agents or aromatic oils as they would damage the product and reduce its life.
- Do not clean the mask in the dishwasher.
- Check carefully all the components to make sure there are no soap deposits or foreign bodies inside the mask.

5 - REGENERATION OF THE MASK FOR REUSE ON DIFFERENT PATIENTS

Notes:

The mask should always be regenerated before being used on a new patient.

- The headgear (e) must be replaced.

The mask can be regenerated for up to 20 cycles by the following processes:

- thermal disinfection at a high level;
- chemical disinfection at a high level;

Detailed mask regeneration instructions can be found on the Internet at www.device.airliquidehealthcare.com or from your retailer.

6 - TROUBLESHOOTING

PROBLEM	POSSIBLE CAUSES	SOLUTIONS
No air comes out of the exhalation orifices.	Flow generator turned off or not working.	<ul style="list-style-type: none"> • Turn on the flow generator; • Adjust the flow generator.
Difficulty breathing.	Patient circuit not inserted correctly.	Connect the patient circuit to the flow generator and the mask.
Noisy mask.	<ul style="list-style-type: none"> • The exhalation orifices are partially clogged; • The mask is assembled incorrectly; • The mask is not positioned correctly. 	<ul style="list-style-type: none"> • Clean the orifices using a soft bristle brush; • Disassemble the mask and then reassemble it as described in the user manual. • Reposition the mask on your face and adjust the headgear as described in the user manual.
Air leaks between the cushion and the face.	<ul style="list-style-type: none"> • The mask is not properly positioned on the face. • The straps of the headgear are not properly adjusted (too tight or too loose). • The size of the mask may not be correct. 	<ul style="list-style-type: none"> • Reposition the mask and adjust the headgear as described in the user manual. • Contact your doctor or the technician to find the correct size


7 - DISPOSAL

The mask components are made of the materials indicated below. They do not contain hazardous substances and can therefore be disposed of with ordinary household waste.

COMPONENT	MATERIAL
Cushion (a)	Silicone + Polycarbonate
Mask body (b)	Polycarbonate
Anti-asphyxia fitting (c)	Polycarbonate + Silicone
Swivelling coupling (d)	Co-polyester
Headgear (e)	Nylon + Polyurethane
Buckle (f)	Nylon

8 - TECHNICAL DATA

PRESSURE CURVE - FLOW

Note: Due to machining differences, the outflow from the exhalation orifices may vary ($\pm 10\%$).

Standards	CE marking in accordance with European Directive 93/42/EEC, Class IIa and subsequent implementation modifications - UNI EN ISO 14971, UNI EN ISO 5356-1, UNI EN ISO 10993-1, UNI EN ISO 17510, UNI EN ISO 17664
Treatment pressure	4 cmH ₂ O – 32 cmH ₂ O
Operating temperature	5°C / 40 °C

ENGLISH

Dead space	69 ml (S) 78 ml (M) 98 ml (L)
Resistance	at 50 l/min: 0.15 cmH ₂ O at 100 l/min: 0.82 cmH ₂ O
A-Weighted Sound Pressure (UNI EN ISO 3744)	12.3 ± 3.2 dB(A)
A-Weighted Sound Power (UNI EN ISO 3744)	20.3 ± 3.2 dB(A)
Overall dimensions maximum (without headgear)	145 mm x 80 mm x 90 mm
Weight (size L)	85 g
Storage temperature	-25°C / +50°C
Storage humidity	0% / 93%
Storage pressure	700 hPa / 1060 hPa

9 - WARRANTY CERTIFICATE

Air Liquide Medical Systems S.r.l. guarantees its product against material or manufacturing defects for 180 days from the date of purchase, provided that the conditions of use in the instructions for use are observed.

If the product is faulty under normal use, Air Liquide Medical Systems S.r.l. will, at its discretion, repair or replace the faulty product or its components.

The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the operating instructions and its intended use. The manufacturer shall be deemed responsible for the safety, reliability and performance of the device if repairs are carried out at the manufacturer or at a laboratory authorized by the manufacturer.

Keep the original casing for the entire warranty period. The warranty in case of product defects can be applied by contacting your dealer.

ENGLISH

For more information on the rights under the warranty, please contact your retailer.

COPYRIGHT

All the information in this manual may not be used for purposes other than the original. This manual is the property of Air Liquide Medical Systems S.r.l. and it may not be reproduced, in whole or in part, without written permission from the company. All rights are reserved.

TECHNICAL CHARACTERISTICS UPDATES

For a continuous improvement of performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l. are regularly subjected to review and change.

Instruction manuals are therefore modified to ensure their constant consistency with the characteristics of devices placed on the market. If the instruction manual accompanying this device is lost, contact the manufacturer to get a copy of the version corresponding to the device by quoting the references shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.