

RESPIREO SOFT FACIAL

Simplicity in comfort



Respireo SOFT Facial



UNIVEO Headgear
Ensures that the mask fits the face optimally


Silicon bumper
Easy regulation of the inclination of the frontal branch

A thin and flexible frame

- For an optimal fit of the cushion on the face thanks to the flexibility of the frontal branch
- Unobstructed field of vision

Cushion SOFT (S, M et L)

- Developed in close partnership with experts
- Balance between flexibility and firmness to ensure optimal and comfortable sealing



A quiet leak
Maximum comfort for the patient and their partner

Ergonomic hooks
Quick assembly and disassembly of the UNIVEO headgear

Ball joint-type
Improves comfort by reducing the tension of the circuit on the mask



Rotating elbow
360° rotation, so as not to disturb the patient in his sleep

Product Reference	Description	Size	Code for order
472548	Respireo SOFT Facial*	S	KM229103
472549	Respireo SOFT Facial*	M	KM229104
472550	Respireo SOFT Facial*	L	KM229105
472551	Respireo SOFT Facial Non Vented *	S	KM229119
472552	Respireo SOFT Facial Non Vented *	M	KM229120
472553	Respireo SOFT Facial Non Vented *	L	KM229121
472633	Respireo SOFT Facial Cushion	S	KM229110
472634	Respireo SOFT Facial Cushion	M	KM229111
472635	Respireo SOFT Facial Cushion	L	KM229112
472471	UNIVEO Headgear**	-	KM227700

► Find all our accompanying videos on Youtube and on our website



Contact

AIR LIQUIDE MEDICAL SYSTEMS
Parc de Haute Technologie
6, rue Georges Besse
92182 Antony Cedex
France

www.device.airliquidehealthcare.com
Mail: dl-frtaema-src-export@airliquide.com

Hotline ALMS: +33 179 51 70 01



Respireo SOFT Facial is a reusable oro-nasal mask with calibrated exhalation orifices (for the Vented version) designed to be used by a single patient at home or by multiple patients in a hospital setting. Respireo SOFT Facial is intended to be used by adult patients (weighing over 30 kg), who have been prescribed non-invasive positive pressure ventilation therapy (NIPPV), e.g. continuous positive airway pressure therapy (CPAP) or Bilevel (positive airway pressure therapy - BiPAP).

Please read the user manual.
Medical device class I CE (**) and IIa CE 0051 (*), Manufactured by Air Liquide Medical Systems S.r.l. Distributed by Air Liquide Medical Systems SA.
Update date: July 2021