

Features and Benefits




Embedded VELCRO®* Tabs enable fast fitting and removal.

F&P Pilairo Q Adjustable Headgear – 400HC326

Adjustable headgear fits a greater range of head sizes and enables sizing adjustment to patient preference.

F&P AirPillow™ Seal – 400HC125

Inflatable AirPillow Seal self-inflates in and around the nose. Can move and flex without compromising the seal.



F&P Pilairo Q Nasal Pillows Mask (no headgear) – 400HC565

Low profile and lightweight Easy Frame for a clear line of sight

Ball-and-Socket Elbow rotates for freedom of movement while reducing drag from the CPAP tube.



F&P Pilairo Q Stretchwise™ Headgear – 400HC322

Minimal Stretchwise™ Headgear automatically adjusts to patients for a quick and easy fit.

F&P Pilairo Q with Adjustable & Stretchwise Heagear – 400421

Mask stabiliser for additional stability

Lightweight flexible tube.

Female Swivel – 400HC224

Swivel with a quiet and smooth rotation.

F&P Pilairo Q Elbow Cover & Diffuser (10 pack) – 400HC229

Q Cover for quieter use¹ and minimal draft.

*VELCRO is a registered trade mark of Velcro Industries B.V.
 1. The Sound Pressure Level of the mask with the Q Cover is 23.6 dBA, with uncertainty 2.5 dBA.